



- Calculating Percent Solutions
 - - Molar Solutions
 - Stock Dilutions

Percent Solutions (expressed as: % = parts per hundred OR % = grams/100 ml)

Weight/Volume When working with solutions for biological applications, distilled or deionized water is usually the solvent of choice. There is the convenient relationship that 1 ml or 1 cc of water weighs 1 gram, so the formula to determine how much solute you need to make an aqueous solution is:

g = (Volume) x (Percent)

Example

You need 550 ml (Volume) of a 7% (Percent) solution of any solute.

Insert the numbers into the formula ...

 $g = (550) \times (0.07)$

g = 38.5 g ... You'll need 38.5 g of solute.

Volume/Volume When 2 or more liquids are mixed together, use this formula:

ml = (Volume) x (%)

Example

You need 75 ml aqueous (aq.) of a 70% ethanol solution (when diluting using 100% ETOH)

Insert the numbers into the formula ...

 $ml = (75) \times (0.7) = 52.5 \, ml$

Measure 52.5 ml of 100 % ethanol

Then bring the volume to 75 ml with dH₂O.

Molar Solutions The formula: grams of solute = (molarity) x (GFM) x (volume in liters)

Example 5

You need to make 250 ml of 0.2 M Cacodylate buffer

Insert the numbers into the formula ...

$$g = (0.2) \times (214) \times (0.250)$$

a = 10.7

Example

You need to make 750 ml of 0.15 M phosphate buffer with a 1:2.5 monobasic: dibasic ratio.

- 1. Determine the volume of each monobasic and dibasic needed.
- 2. Take the final volume (in this case 750 ml) and divide it by the total parts of the ratio:

3. The volume for the monobasic part is:

$$1 \times 214.3 \text{ ml} = 214.3 \quad g = (0.15) \times (138) \times (0.2143) \quad g = 4.43$$

4. The volume for the dibasic part is:

$$2.5 \times 214.3 \text{ ml} = 535.75$$
 $q = 0.15) \times (268) \times (0.5357)$ $q = 21.54$

Summation: Add 4.43 grams of monobasic and 21. 54 grams of dibasic for a final volume of 750 ml

Stock Dilutions The formula is C₁ x V₁ = C₂ x V₂ ...

C₁ is the concentration of the stock solution

 $\mathbf{V_1}$ is the volume of the stock solution, usually the unknown part of the equation

C, is the concentration of the working solution

V₂ is the volume of the working solution

Example You pr You need to make a working solution of 250 ml of 2.5% glutaraldehyde using a 50% stock solution 50% stock solution.

Insert the numbers into the formula ...

$$C_1$$
 (50%) x V_1 (?) = C_2 (2.5%) x V_2 (250)

... solve for V₁ ...

 $V_1 = 625/50 = 12.5$ ml of 50% into a volume of 250 ml

Example

Using the above example, you have a 10 ml ampule of 50% glutaraldehyde and you want to use it all.

Insert the numbers into the formula ... 9

$$C_1 (50\%) \times V_1 (10 \text{ ml}) = C_2 (2.5\%) \times V_2 (? \text{ ml})$$

... solve for V₂ ...

 $V_{2} = 500/2.5 = 200 \text{ ml}$ is the final volume

Example

This example uses a Molar solution instead of a Percent solution: You need a working solution – 350 ml of 0.125 M buffer using a 2 M stock solution.

Insert the numbers into the formula ...

ctron Microscopy Sciences

$$C_1 (2 \text{ M}) \times V_1 (?) = C_2 (0.125 \text{M}) \times V_2 (350)$$

solve for V₁ ...

 $V_1 = 43.75/2 = 21.875$ ml of 2 M stock into a volume of 350 ml